



HIV type 1 p17 (IIIB) [FITC] (DAGH006)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate FITC

Format Frozen Liquid

Preservative None

BACKGROUND

Keywords HIV type 1 p17 ; HIV-1 gp17; HIV1 gp17; Envelope surface glycoprotein gp17; Glycoprotein 17; gp17; gp17 glycoprotein; Human Immunodeficiency Virus 1; SU; Surface protein; Retroviridae; Lentivirus; human immunodeficiency virus