



Mouse Anti-Human PCT monoclonal antibody, clone IN089 (DMABP-L38)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified recombinant human Procalcitonin
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	IN089
Purification	Protein A Purified
Conjugate	Unconjugated
Applications	LFIA. We recommend the following for sandwich ELISA (Capture - Detection): DMABP-L38 - DMABP-L39 DMABP-L38 - CABT-L6523Z
Epitope	N-term PCT
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.0
Preservative	0.09% Sodium Azide

Storage Store at 4°C for short term. Aliquot and store at -20°C for long term. Avoid freeze/thaw cycles.

Ship Wet ice

BACKGROUND

Introduction Procalcitonin (PCT) is a peptide precursor of the hormone calcitonin, the latter being involved with calcium homeostasis. It arises once preprocalcitonin is cleaved by endopeptidase. It was first identified by Leonard J. Deftos and Bernard A. Roos in the 1970s. It is composed of 116 amino acids and is produced by parafollicular cells (C cells) of the thyroid and by the neuroendocrine cells of the lung and the intestine.

Keywords Procalcitonin; PCT; PCT antigen; PCT protein
