



Goat Anti-Echis ocellatus Venom Polyclonal Antibody (DPAB-JX2371342)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Immunogen | Echis ocellatus Venom |
| Isotype | IgG |
| Source/Host | Goat |
| Species Reactivity | N/A |
| Purification | Unpurified |
| Conjugate | Unconjugated |
| Applications | IA Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Format | Unpurified, antiserum |
| Size | 1 mL |
| Preservative | None |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Ship | Wet ice |